

# RAPID RELIEF WRIST

## THANK YOU FOR YOUR PURCHASE OF COPPER FIT®

Copper Fit's Rapid Relief Wrist features a "one size fits most" brace that holds a Hot / Cold Gel pack in an interior pouch. Hot / Cold Therapy reduces muscle/joint pain and inflammation. Rapid Relief Braces/Wraps are designed to conform to your body, letting you move about freely while administering Hot / Cold Therapy. Read instructions thoroughly before using Copper Fit Rapid Relief Braces/Wraps and Hot /Cold Gel packs.

Visit [www.clikit.ca](http://www.clikit.ca) for the full collection

### CAUTION:

Wash product prior to use. Product should be snug, but not so tight that it restricts circulation. Tight fitting garments can cause heat rash or other skin conditions, especially when hot and sweaty. Discontinue use if rash develops. If these symptoms occur, discontinue use and consult a physician. Do not place Rapid Relief Wrap in freezer, microwave or hot water; only gel pack.

### NEOPRENE:

Product contains synthetic rubber Neoprene (polychloroprene). Do not use if you have an allergy to Neoprene.

#### RAPID RELIEF WRAP CARE AND MAINTENANCE:

HAND WASH IN COLD WATER • AIR DRY • DO NOT BLEACH • DO NOT IRON  
• DO NOT DRY CLEAN • DO NOT TUMBLE DRY

The gel pack should be used in 20 minute intervals. Applying the gel pack for longer than the suggested time is not recommended.

### TO USE COLD:

Place the gel pack flat in the freezer for at least 2 hours prior to use. Moisture in the gel pack will cause frosting. Gel will remain cold even after frost is no longer visible.

### TO USE HEATED:

#### HOT WATER METHOD:

CAUTION: NEVER TOUCH THE GEL PACK WITH BARE HANDS WHEN SUBMERGING/RETRIEVING FROM HOT WATER. Use an appropriate tool such as a wooden spoon. Carefully submerge the gel pack into a bowl of hot water for 20-30 seconds. Using the tool, carefully remove the gel pack from the bowl and place onto a towel. If you find the gel pack is not at your desired temperature, carefully submerge the gel pack for an additional 10 seconds at a time until desired temperature is reached.

#### MICROWAVE METHOD:

Wrap your gel pack with a damp cloth or paper towel. Lay flat in the microwave and heat as follows:

|                         |            |            |            |
|-------------------------|------------|------------|------------|
| Suggested Heating Times | 800 Watts  | 1000 Watts | 1200 Watts |
|                         | 40 Seconds | 30 Seconds | 20 Seconds |

If you find the gel pack is not at the desired temperature, flip over the gel pack and microwave an additional 10 seconds until desired temperature is reached. DO NOT exceed 2 (two) 10 second additional cycles. Heating the gel pack beyond the advised time is not recommended.

WARNING! Gel pack and towel will be hot when removing from the microwave.

WARNING! At no time should the gel pack be in use if it exceeds 49°C (120°F)

CAUTION: This product can cause burns or frostbite. DO NOT place gel pack directly on skin. DO NOT use on infants, children or those unable to remove the product. Consult a physician prior to use if you have diabetes, poor circulation, rheumatoid arthritis or other chronic conditions. DO NOT use if you are pregnant. DO NOT use on sensitive skin or insensitve (numb) skin or on an open wound. DO NOT wear while sleeping. DO NOT swallow gel. DO NOT open gel pack and discontinue use if gel pack is punctured, leaking or damaged in any way. DO NOT place any heavy objects or extra pressure on top of the gel pack. If irritation or discomfort occurs, remove product immediately. Stop use and consult a physician if pain worsens or persists for more than 3 days. KEEP OUT OF REACH OF CHILDREN AND PETS

CARE: Gel pack surface can be cleaned with mild soap, detergent or alcohol.

### LIMITATION OF LIABILITY

Liability limited to the purchase price of this product. Supertek Canada Inc. shall not be liable for any incidental or consequential damages or any implied warranty on this product.

#### Distributed by:

SUPERTEK CANADA Inc., Mont-Royal QC H4T 1X2

[WWW.CLIKIT.CA](http://WWW.CLIKIT.CA)


Customer Service: 1-800-304-8354 • For more information, visit [clikit.ca](http://clikit.ca)

# RAPID RELIEF WRIST


The Rapid Relief Wrist is designed to be worn on either the right or left wrist.  
 NOTE: The Rapid Relief Wrist Brace can be used with or without the gel pack.

## HOW TO USE

1. Prepare the Hot/Cold gel pack according to directions printed on the gel pack.
2. Undo all straps on the wrist brace.
3. Place wrist brace on flat surface with inside facing up (Image A).
4. To insert Hot/Cold gel pack, undo the fasteners found on the inside bottom part of the wrist brace (Image A).
5. Slide the gel pack into the interior pouch with the long tips pointing upward (Images B & C).
6. Secure fasteners together to close the interior pouch opening.



7. Use the wrist brace on either hand by following the right and left labels below the thumb holes.
8. Place your desired wrist in the middle of the wrist brace and place your thumb in the corresponding thumb hole (Image D).
9. Fold over each side of the wrist brace to feed the straps through (Image E).
10. Feed bottom strap through the buckle and around to secure (Image F).
11. Feed top strap through buckle and secure to desired fit (Image G).
12. Adjust each strap to customize your fit.



# RAPID RELIEF WRIST

## THANK YOU FOR YOUR PURCHASE OF COPPER FIT®

Copper Fit's Rapid Relief Wrist features a "one size fits most" brace that holds a Hot / Cold Gel pack in an interior pouch. Hot / Cold Therapy reduces muscle/joint pain and inflammation. Rapid Relief Braces/Wraps are designed to conform to your body, letting you move about freely while administering Hot / Cold Therapy. Read instructions thoroughly before using Copper Fit Rapid Relief Braces/Wraps and Hot /Cold Gel packs.

Visit [www.clikit.ca](http://www.clikit.ca) for the full collection

### CAUTION:

Wash product prior to use. Product should be snug, but not so tight that it restricts circulation. Tight fitting garments can cause heat rash or other skin conditions, especially when hot and sweaty. Discontinue use if rash develops. If these symptoms occur, discontinue use and consult a physician. Do not place Rapid Relief Wrap in freezer, microwave or hot water; only gel pack.

### NEOPRENE:

Product contains synthetic rubber Neoprene (polychloroprene). Do not use if you have an allergy to Neoprene.

#### RAPID RELIEF WRAP CARE AND MAINTENANCE:

HAND WASH IN COLD WATER • AIR DRY • DO NOT BLEACH • DO NOT IRON  
• DO NOT DRY CLEAN • DO NOT TUMBLE DRY

The gel pack should be used in 20 minute intervals. Applying the gel pack for longer than the suggested time is not recommended.

### TO USE COLD:

Place the gel pack flat in the freezer for at least 2 hours prior to use. Moisture in the gel pack will cause frosting. Gel will remain cold even after frost is no longer visible.

### TO USE HEATED:

#### HOT WATER METHOD:

CAUTION: NEVER TOUCH THE GEL PACK WITH BARE HANDS WHEN SUBMERGING/RETRIEVING FROM HOT WATER. Use an appropriate tool such as a wooden spoon. Carefully submerge the gel pack into a bowl of hot water for 20-30 seconds. Using the tool, carefully remove the gel pack from the bowl and place onto a towel. If you find the gel pack is not at your desired temperature, carefully submerge the gel pack for an additional 10 seconds at a time until desired temperature is reached.

#### MICROWAVE METHOD:

Wrap your gel pack with a damp cloth or paper towel. Lay flat in the microwave and heat as follows:

|                         |            |            |            |
|-------------------------|------------|------------|------------|
| Suggested Heating Times | 800 Watts  | 1000 Watts | 1200 Watts |
|                         | 40 Seconds | 30 Seconds | 20 Seconds |

If you find the gel pack is not at the desired temperature, flip over the gel pack and microwave an additional 10 seconds until desired temperature is reached. DO NOT exceed 2 (two) 10 second additional cycles. Heating the gel pack beyond the advised time is not recommended.

WARNING! Gel pack and towel will be hot when removing from the microwave.

WARNING! At no time should the gel pack be in use if it exceeds 49°C (120°F)

CAUTION: This product can cause burns or frostbite. DO NOT place gel pack directly on skin. DO NOT use on infants, children or those unable to remove the product. Consult a physician prior to use if you have diabetes, poor circulation, rheumatoid arthritis or other chronic conditions. DO NOT use if you are pregnant. DO NOT use on sensitive skin or insensitive (numb) skin or on an open wound. DO NOT wear while sleeping. DO NOT swallow gel. DO NOT open gel pack and discontinue use if gel pack is punctured, leaking or damaged in any way. DO NOT place any heavy objects or extra pressure on top of the gel pack. If irritation or discomfort occurs, remove product immediately. Stop use and consult a physician if pain worsens or persists for more than 3 days. KEEP OUT OF REACH OF CHILDREN AND PETS

CARE: Gel pack surface can be cleaned with mild soap, detergent or alcohol.

### LIMITATION OF LIABILITY

Liability limited to the purchase price of this product. Supertek Canada Inc. shall not be liable for any incidental or consequential damages or any implied warranty on this product.

#### Distributed by:

SUPERTEK CANADA Inc., Mont-Royal QC H4T 1X2

[WWW.CLIKIT.CA](http://WWW.CLIKIT.CA)

Customer Service: 1-800-304-8354 • For more information, visit [clikit.ca](http://clikit.ca)


# RAPID RELIEF WRIST

The Rapid Relief Wrist is designed to be worn on either the right or left wrist.

NOTE: The Rapid Relief Wrist Brace can be used with or without the gel pack.

## HOW TO USE

1. Prepare the Hot/Cold gel pack according to directions printed on the gel pack.
2. Undo all straps on the wrist brace.
3. Place wrist brace on flat surface with inside facing up (Image A).
4. To insert Hot/Cold gel pack, undo the fasteners found on the inside bottom part of the wrist brace (Image A).
5. Slide the gel pack into the interior pouch with the long tips pointing upward (Images B & C).
6. Secure fasteners together to close the interior pouch opening.



7. Use the wrist brace on either hand by following the right and left labels below the thumb holes.
8. Place your desired wrist in the middle of the wrist brace and place your thumb in the corresponding thumb hole (Image D).
9. Fold over each side of the wrist brace to feed the straps through (Image E).
10. Feed bottom strap through the buckle and around to secure (Image F).
11. Feed top strap through buckle and secure to desired fit (Image G).
12. Adjust each strap to customize your fit.

